

## Middleby Cooking Systems Group



1400 Toastmaster Drive • Elgin, IL 60120-9272 U. S. A. • Telephone: (847)741-3300 • Service and Parts Fax: (847)741-4406

## Middleby Marshall® Service Bulletin

### SUBJECT

Availability of center lift tube for stacking PS555 oven cavities

### GENERAL

A new center lift tube, P/N 42966, is now available. This tube can be used with the existing PS570 stacking lift kit (P/N 30580) in order to stack PS555 oven cavities. The two lift tubes in the kit (P/Ns 30557) will have to be modified (as described below) to allow the new center tube to be used. After this modification, the stacking lift kit can be used to stack both PS555 and PS570 oven cavities.

### APPLICATION

In order to be able to stack PS555 ovens, first order one center lift tube, P/N 42966, from the Middleby parts department. If you do not have the PS570 stacking lift kit, it can also be ordered at this time (P/N 30580).

To adapt the stacking lift kit for use with the new center tube, refer to the drawing on Page 2 of this Bulletin. Cut 7-1/2" (191mm) from the inner end of each lift tube (P/N 30557). Note that the inner end is the end without the eight 5/16" (8mm) holes.

After this modification, the lift kit can be used to stack PS555 oven cavities by using the new center tube. To stack PS570 oven cavities, use the kit's original center lift tube, P/N 30596.

### COMPATIBILITY

After modifying the kit's two lift tubes, and ordering the new center lift tube, the PS570 Stacking Lift Kit is compatible with BOTH PS570 and PS555 ovens.

### AVAILABILITY

The PS555 center lift tube, P/N 42966, is now available and can be ordered normally. This part is offered to our authorized service agencies at a cost of \$15.36, plus shipping and handling.

